

## Sartolab® RF|BT

### Vacuum Filtration units



Sartolab® RF



Sartolab® RF 50

- Sartolab RF와 Sartolab BT는 Cell Culture시, 최대의 효율 보장
- 빠른 유량과 적은 Binding을 유도하는 PES Membrane 장착
- Tube Adapter와 Lid 포함
- 여러 종류의 Volume Size를 가짐

### Ordering Information for Sartolab® RF Vacuum Filtration Units

Cat. No.	Description	Membrane Type	Pore Size	Funnel Volume(ml)	Bottle Unit
180E01-2	Sartolab® RF 50	Asymmetric PES	0.22μm	50	24
180F01-2	Sartolab® RF 50	PES	0.45μm	50	24
180E01-E8	Sartolab® RF 50	Asymmetric PES	0.22μm	50	96*
180F01-E8	Sartolab® RF 50	PES	0.45μm	50	96*
180E02-E	Sartolab® RF 150	Asymmetric PES	0.22μm	150	12
180F02-E	Sartolab® RF 150	PES	0.45μm	150	12
180D03-E	Sartolab® RF 250	PES	0.1μm	250	12
180E03-E	Sartolab® RF 250	Asymmetric PES	0.22μm	250	12
180F03-E	Sartolab® RF 250	PES	0.45μm	250	12
180E04-E	Sartolab® RF 500	Asymmetric PES	0.22μm	500	12
180F04-E	Sartolab® RF 500	PES	0.45μm	500	12
180D05-E	Sartolab® RF 1,000	PES	0.1μm	1,000	12
180E05-E	Sartolab® RF 1,000	Asymmetric PES	0.22μm	1,000	12
180F05-E	Sartolab® RF 1,000	PES	0.45μm	1,000	12

\* 96 pack: non-sterile, without tubing connector



Sartolab® BT

### Ordering Information for Sartolab® BT Vacuum Filtration Units

Cat. No.	Description	Membrane Type	Pore Size	Funnel Volume(ml)	Bottle Unit
180E12-E	Sartolab® BT 150	Asymmetric PES	0.22μm	150	12
180E13-E	Sartolab® BT 250	Asymmetric PES	0.22μm	250	12
180E14-E	Sartolab® BT 500	Asymmetric PES	0.22μm	500	12
180E15-E	Sartolab® BT 1,000	Asymmetric PES	0.22μm	1,000	12
180F15-E	Sartolab® BT 1,000	PES	0.45μm	1,000	12



Sartolab® bottle

### Ordering Information for Sartolab® bottle Units

#### Sartolab® Bottle

Delivered sterile, for filtration and storage

Cat. No.	Description	Volume(ml)	Packing Unit
180-22-E	Sartolab® bottle 150ml	150	12
180-23-E	Sartolab® bottle 250ml	250	12
180-24-E	Sartolab® bottle 500ml	500	12
180-25-E	Sartolab® bottle 1,000ml	1,000	12